

- ○ QC Wireless SDK Instructions for Use
  - ■ ■ update note:
  - 1.introduce
    - ■ 1.1The role of the SDK
  - 2. API description
  - 2.0 Access Conditions
  - 2.1 Permissions required by the SDK
  - 2.2 Access conditions
  - 2.3.1API
    - Scanning Devices
    - device connection:BleOperateManager.getInstance()
  - 2.3.2 Feature list:
    - Synchronize time, get the list of functions supported by the device
    - bracelet battery
    - Continuous heart rate, blood oxygen, blood pressure switch
    - Set watch sports goals
    - Find equipment
    - Ring photo control
    - Set the Ring to factory reset
    - firmware version number, hardware version number
  - 2.3.3 Data synchronization:
    - Synchronize steps, distance, kcal for the day
    - Synchronized step data details
    - Sync Sleep Data Details
    - Synchronize the details of new sleep data and return according to SetTimeRsp
    - Sync heart rate data
    - Synchronized blood pressure function
    - Synchronized blood oxygen function
    - Synchronized pressure function
  - 2.3.6 OTA upgrade function:
  - 2.3.7Manual measurement

- 2.3.9 Changes in Ring measurement data are proactively reported to the APP
- 2.3.10 APP opens exercise type

## QC Wireless SDK Instructions for Use

---

1. Author: James
2. Shenzhen QC.wireless Technology Co., Ltd.
3. Version: 1.0.0

### update note:

1. (2021/07/06) scan, connect, measure commands
2. (2021/07/20) Add setting command
3. (2021/07/21) Increase step count, heart rate, sleep data sync
4. (2022/02/28) Add new sleep algorithm
5. (2023/03/10) Add message switch synchronization, body temperature, user information, watch manual measurement of heart rate, blood pressure results
6. (2023/03/16) Add bt connection and contact person
7. (2024/02/23) Add manual pressure, APP movement, pressure synchronization, pressure setting

## 1.introduce

### 1.1The role of the SDK

Provide partner companies with the Android Bluetooth SDK for use with Green Orange wireless devices that provide basic and advanced functionality for a major watch or other device. This document is intended to explain the usage context, functionality, etc. of the API. Intended Audience and Reading Recommendations The intended audience and reader recommendations in this article are shown in Annex 1.

| Reader                         | Role                                                                                      |
|--------------------------------|-------------------------------------------------------------------------------------------|
| Software Architecture Engineer | Architecture Analysis and Technical Guidance                                              |
| Android development engineer   | Have a certain android development ability, understand Ble related development technology |

## 2. API description

### 2.0 Access Conditions

Android 5.0 or above, Bluetooth 4.0 or above.

### 2.1 Permissions required by the SDK

```

//network permissions
<uses-permission android:name="android.permission.INTERNET" />
//Bluetooth related permissions
<uses-permission android:name="android.permission.BLUETOOTH" />
<uses-permission
android:name="android.permission.BLUETOOTH_ADMIN" />
//Storage related permissions
<uses-permission
android:name="android.permission.WRITE_EXTERNAL_STORAGE" />
<uses-permission
android:name="android.permission.READ_EXTERNAL_STORAGE" />
//Location permission
<uses-permission
android:name="android.permission.ACCESS_COARSE_LOCATION" />
<uses-permission
android:name="android.permission.ACCESS_BACKGROUND_LOCATION" />
<uses-permission
android:name="android.permission.ACCESS_FINE_LOCATION" />

If it is android 12 or higher system
<uses-permission
android:name="android.permission.BLUETOOTH_CONNECT" />
<uses-permission
android:name="android.permission.BLUETOOTH_SCAN" />
<uses-permission
android:name="android.permission.BLUETOOTH_ADVERTISE" />

```

## 2.2 Access conditions

- Green Orange Wireless Wearables
- Green Orange Wireless SDK and documentation

### 2.3.1API

#### Scanning Devices

```

        //start scan
        BleScannerHelper.getInstance().scanDevice(final Context
context, UUID mUuid, final ScanWrapperCallback scanCallBack);
        //stop scan
        BleScannerHelper.getInstance().stopScan(Context context)
        //Specified device scan
        BleScannerHelper.getInstance().scanTheDevice(final Context
context, final String macAddress, final OnTheScanResult scanResult)

```

## device connection:BleOperateManager.getInstance()

```

        //direct connection

        BleOperateManager.getInstance().connectDirectly(smartWatch.deviceAd
dress)
        //scan connections

        BleOperateManager.getInstance().connectWithScan(smartWatch.deviceAd
dress)
        //disconnect
        BleOperateManager.getInstance().unBindDevice()
        //reconnect
        BleOperateManager.getInstance().setNeedConnect(boolean
needConnect)
        //Called when bluetooth is turned off
        BleOperateManager.getInstance().setBluetoothTurnOff(false)
        BleOperateManager.getInstance().disconnect()
        //Turn on the system bluetooth monitor
        BleOperateManager.getInstance().setBluetoothTurnOff(true)

```

## 2.3.2 Feature list:

Synchronize time, get the list of functions supported by the device

```

        //set time
        CommandHandle.getInstance().executeReqCmd(SetTimeReq(0),
        ICommandResponse<SetTimeRsp>() {})
        //Callback Description
        public class SetTimeRsp extends BaseRspCmd {
        //support body temperature
        public boolean mSupportTemperature;
        //watch face
        public boolean mSupportPlate;
        //Support the menstrual cycle
        public boolean mSupportMenstruation;
        //Support custom watch faces
        public boolean mSupportCustomWallpaper;
        //Support blood oxygen
        public boolean mSupportBloodOxygen;
        //blood pressure support
        public boolean mSupportBloodPressure;
        //Support fatigue
        public boolean mSupportFeature;
        //Support one-key detection
        public boolean mSupportOneKeyCheck;
        //weather support
        public boolean mSupportWeather;
        //Support for new sleep protocol
        public boolean mNewSleepProtocol;
        //Supports up to 6 or 3 dials
        public int mMaxWatchFace;
    }

```

**bracelet battery**

```

        CommandHandle.getInstance().executeReqCmd(new
SimpleKeyReq(Constants.CMD_GET_DEVICE_ELECTRICITY_VALUE), new
ICommandResponse<BatteryRsp>() {

            @Override
            public void onDataResponse(BatteryRsp resultEntity) {
                if (resultEntity.getStatus() ==
BaseRspCmd.RESULT_OK) {
                    //Get battery successfully
                }
            }

        });

//Callback Description
public class BatteryRsp extends BaseRspCmd {
//battery [0-100]
private int batteryValue;
}

```

## Continuous heart rate, blood oxygen, blood pressure switch

```

//Read continuous heart rate settings
CommandHandle.getInstance()
    .executeReqCmd(HeartRateSettingReq.getReadInstance(),
        ICommandResponse<HeartRateSettingRsp> {
            //it.isEnabled        switch
            //it.heartInterval    heart rate measurement interval

        })

//Read Continuous SpO2 settings
CommandHandle.getInstance()
    .executeReqCmd(BloodOxygenSettingReq.getReadInstance(),
        ICommandResponse<BloodOxygenSettingRsp> {

        })

//Read continuous blood pressure settings
CommandHandle.getInstance()
    .executeReqCmd(BpSettingReq.getReadInstance(),
        ICommandResponse<BpSettingRsp> {

        })

//Read pressure setting
CommandHandle.getInstance()

```

```

        .executeReqCmd(PressureSettingReq.getReadInstance(),
            ICommandResponse<PressureSettingRsp>() {
                switch: it.isEnabled
            })

        //Write continuous heart rate switch isEnabled: true on,
        false: off  hrInterval: 10, 15, 20, 30, 60

        CommandHandle.getInstance().executeReqCmd(
            HeartRateSettingReq.getWriteInstance(true, hrInterval),
            ICommandResponse<HeartRateSettingRsp> {

                })
        //Write continuous blood oxygen switch isEnabled: true on,
        false: off
        CommandHandle.getInstance().executeReqCmd(
            BloodOxygenSettingReq.getWriteInstance(boolean isEnabled),
            ICommandResponse<BloodOxygenSettingRsp> {

                })
        //write blood pressure switch

        CommandHandle.getInstance().executeReqCmd(BpSettingReq.getWriteInst
            ance(boolean isEnabled, StartEndTimeEntity startEndTimeEntity, int
            multiple),  ICommandResponse<BpSettingRsp> {

                })
        BpSettingRsp,Parameter Description
        isEnabled: true on false off
        StartEndTimeEntity The parameter description is the same as
        above
        multiple  default 60

        //Write pressure setting switch,switch isEnabled: true on,
        false: off
        CommandHandle.getInstance().executeReqCmd(
            PressureSettingReq.getWriteInstance(isEnabled),
            ICommandResponse<PressureSettingRsp> {

                })
    }
}

```

## Set watch sports goals



```
CommandHandle.getInstance().executeReqCmd(  
    TargetSettingReq.getWriteInstance(  
        final int step, final int calorie, final int  
distance, final int sportMinute, final int sleepMinute  
    ), null  
)
```

#### Parameter Description

step: step target

calorie: Calorie target target, unit card, write kcal to a\*1000

distance: Distance to target, in meters

sportMinute: Exercise duration target, in minutes

sleepMinute: Sleep duration target, in minutes

## Find equipment

```
CommandHandle.getInstance().executeReqCmd(FindDeviceReq(), null)
```

## Ring photo control

```

//The bracelet enters the camera interface

CommandHandle.getInstance().executeReqCmd(CameraReq(CameraReq.ACTION_INT0_CAMARA_UI), null)
//The wristband is controlled by the bright screen on the camera interface. The APP will send the bright screen command to keep the watch bright. It is recommended to send it every 2 seconds.
CommandHandle.getInstance().executeReqCmd(

CameraReq(CameraReq.ACTION_KEEP_SCREEN_ON),
        null
    )
//Bracelet click to take a photo event monitoring

BleOperateManager.getInstance().addNotifyListener(Constants.CMD_TAKING_PICTURE,new ICommandResponse<CameraNotifyRsp>(){

    @Override
    public void onDataResponse(CameraNotifyRsp
resultEntity) {

        }
    });
resultEntity.getAction()
Parameter Description
//The watch exits the camera interface
CameraNotifyRsp.ACTION_FINISH
//The watch clicked on the photo event
CameraNotifyRsp.ACTION_TAKE_PHOTO

CommandHandle.getInstance().executeReqCmd(CameraReq(CameraReq.ACTION_FINISH), null)

```

## Set the Ring to factory reset

```

CommandHandle.getInstance().executeReqCmd(RestoreKeyReq(Constants.CMD_RE_STORE),null)

```

- message push

```
//The watch message push switch is fully turned on, and the APP  
should be actively opened  
    CommandHandle.getInstance().executeReqCmd(  
        SetANCSReq(), null  
    )  
//Send message push to watch  
    MessPushUtil.pushMsg(type,message:String)
```

PushMsgUintReq parameter description  
type:

- 0x00: Call reminder 0x01: SMS reminder 0x02: QQ reminder
- 0x03: WeChat reminder,
- 0x04: incoming call to answer or hang up 0x05: Facebook  
message reminder 0x06: WhatsApp message reminder
- 0x07: Twitter message reminder 0x08: Skype message reminder
- 0x09: Line message reminder 0x0a: LinkedIn
- 0x0b: Instagram 0x0c: TIM message 0x0d: Snapchat
- 0x0e: others other types of notifications

**firmware version number, hardware version number**

```

        //hardware information

CommandHandle.getInstance().execReadCmd(CommandHandle.getInstance()
    .getReadHwRequest());
        //firmware information

CommandHandle.getInstance().execReadCmd(CommandHandle.getInstance()
    .getReadFmRequest());

    Receiving implements this QCBluetoothCallbackCloneReceiver
    refer to demo MyBluetoothReceiver
    Judging UUID in the callback onCharacteristicRead

    override fun onCharacteristicRead(uuid: String?, data:
    ByteArray?) {
        if (uuid != null && data != null) {
            val version = String(data, Charsets.UTF_8)
            when(uuid){
                Constants.CHAR_FIRMWARE_REVISION.toString() -> {
                    //Firmware version number version
                }
                Constants.CHAR_HW_REVISION.toString() -> {
                    //hardware version number version number
version
                }
            }
        }
    }
}

```

### 2.3.3 Data synchronization:

Synchronize steps, distance, kcal for the day

```
CommandHandle.getInstance().executeReqCmd(  
    SimpleKeyReq(Constants.CMD_GET_STEP_TODAY),  
    ICommandResponse<TodaySportDataRsp> {})
```

TodaySportDataRsp parameter description

```
// days ago  
private int daysAgo;  
// date: year  
private int year;  
// date: month  
private int month;  
// date: day  
private int day;  
// total steps  
private int totalSteps;  
// running steps/aerobic steps  
private int runningSteps;  
// calorie value  
private int calorie;  
// walking distance  
private int walkDistance;  
// Movement time, in seconds  
private int sportDuration;  
// sleep time in seconds  
private int sleepDuration;
```

## Synchronized step data details

```

    dayOffset 0: Today 1: Yesterday 2: The day before yesterday,
supports synchronization for up to 7 days
    CommandHandle.getInstance().executeReqCmd(
        ReadDetailSportDataReq(dayOffset, 0, 95),
        ICommandResponse<ReadDetailSportDataRsp> {

            })
BleStepDetails parameter description
    //year
    private int year;
    //moon
    private int month;
    //day
    private int day;
    //15 minutes a point, the total number of points in a day is
96 points, [0, 95], used to calculate the details of the number of
steps per hour
    private int timeIndex=0;
    // calorie unit card
    private int calorie=0;
    //Step count
    private int walkSteps=0;
    //distance in meters
    private int distance=0;
    // keep for now
    private int runSteps=0;

```

## Sync Sleep Data Details

```

    //deviceAddress Device mac address
    //dayOffset 0: Today 1: Yesterday 2: The day before yesterday,
    supports synchronization for up to 7 days
    //ISleepCallback callback

SleepAnalyzerUtils.getInstance().syncSleepReturnSleepDisplay(device
Address,dayOffset, ISleepCallback {
    //callback: SleepDisplay
})

SleepDisplay parameter description
    // total sleep time
    private int totalSleepDuration;
    // total time of deep sleep
    private int deepDuration;
    // total time of light sleep
    private int shallowDuration;
    // Go to sleep timestamp in seconds
    private int sleepTime;
    // wake up timestamp in seconds
    private int wakeTime;
    // A set of sleep data
    private List<SleepDataBean> list;
    private String address;

SleepDataBean parameter description
    sleepStart start timestamp of a sleep in seconds
    sleepEnd The end timestamp of a sleep in seconds
    type sleep type 2: light sleep 1: deep sleep 3: awake

```

**Synchronize the details of new sleep data and return according to SetTimeRsp**

```

    //offset 0 today 1 yesterday
    public void syncSleepList(int offset, final
ILargeDataSleepResponse response)
    SleepNewProtoResp    parameter description
    //sleep start time
    private int st;
    //sleep end time
    private int et;
    //sleep collection
    private List<DetailBean> list;
    DetailBean parameter description
    // duration of a sleep type
    private int d;
    //type of sleep 2: light sleep 3: deep sleep 5: awake
    private int t;

```

## Sync heart rate data

```

nowTime current time zone * 3600 + unix second value of current
time
    Sync yesterday: nowTime-(24*3600)*1,
    Sync the day before yesterday: nowTime-(24*3600)*2
    Data can be synchronized for up to three days
    val time = (getTimeZone() * 3600).toInt()
    val nowTime = date.unixTimestamp + time

    CommandHandle.getInstance().executeReqCmd(
        ReadHeartRateReq(nowTime),
        ICommandResponse<ReadHeartRateRsp> {

        })

ReadHeartRateRsp parameter description
    //nothing yet
    private int size = 0;
    //nothing yet
    private int index = 0;
    // unix second value of heart rate data
    private int mUtcTime;
    //The heart rate data array is one point every 5 minutes, the
    data subscript *5 is equal to the number of minutes of the day
    private byte[] mHeartRateArray;
    private boolean endFlag = false;

```



## Synchronized blood pressure function

```
//Synchronized automatic blood pressure, measured once an hour
CommandHandle.getInstance()

.executeReqCmd(SimpleKeyReq(Constants.CMD_BP_TIMING_MONITOR_DATA),
ICommandResponse<BpDataRsp> {})
    BpDataEntity parameter description
    //year
    private int year;
    //month
    private int month;
    //day
    private int day;
    private int timeDelay;
    private ArrayList<BpValue> bpValues;

    BpValue parameter description
    //The minute of the day, usually the whole hour
    int timeMinute;
    //measured heart rate value
    int value;

    Get the blood pressure value calculated from the measured value
    //The heart rate value returned by the hr callback, age is the
    age of the user
    val sbp= CalcBloodPressureByHeart.cal_sbp(hr, age) (systolic
    blood pressure)
    //sbp heart rate calculated value
    val dbp=CalcBloodPressureByHeart.cal_dbp(sbp) (diastolic
    pressure)

    //Confirm blood pressure synchronization, call after receiving
    the callback, the watch will delete the records that have been
    synchronized

CommandHandle.getInstance().executeReqCmd(BpReadConformReq(true),nu
ll)

    //synchronize manual blood pressure
    CommandHandle.getInstance()
        .executeReqCmd(ReadPressureReq(0),
ICommandResponse<ReadBlePressureRsp> {})

    ReadBlePressureRsp.getValueList() parameter description
    BpPressure parameter description
    //time seconds value
    public long time;
```

```

public long time;
//(Diastolic pressure)

public int dbp;
//(systolic blood pressure)
public int sbp;

```

## Synchronized blood oxygen function

```

    LargeDataHandler.getInstance().syncBloodOxygenWithCallback(new
IBloodOxygenCallback() {
        @Override
        public void readBloodOxygen(List<BloodOxygenEntity>
data) {

            }
        });

    BloodOxygenEntity parameter description
        //data date
        private String dateStr;
        //Data minimum value array, one data per hour, a total of
24
        private List<Integer> minArray;
        //Data maximum value array, one data per hour, a total of
24
        private List<Integer> maxArray;
        //Data value at 0:00 on a certain day
        private long unix_time;

```

## Synchronized pressure function

```

        //offset 0-6
        //Synchronization pressure 7 days 0 only synchronizes today 1
        yesterday 2 synchronizes the day before yesterday .... supports up
        to 7 days
        CommandHandle.getInstance()
            .executeReqCmd(PressureReq(offset),
                ICommandResponse<PressureRsp> {

            })

        //PressureRsp Parameter Description
        //Pressure data, the return value divided by 10 is the value
        displayed by the APP, and one data is generated every half hour.
        private byte[] pressureArray;
        //Stress test interval
        private int range=30;
        //Pressure data date
        private DateUtil today;

```

### 2.3.6 OTA upgrade function:

```

//dfu upgrade instance
val fuHandle= DfuHandle.getInstance()
//initialize callback
dfuHandle.initCallback()
//DFU file verification, path firmware file path
if (dfuHandle.checkFile(path)) {
    dfuHandle.start(dfuOpResult)
}
//dfuOpResult callback description
private val dfuOpResult: DfuHandle.IOpResult = object :
DfuHandle.IOpResult {
    override fun onActionResult(type: Int, errCode: Int) {
        if (errCode == DfuHandle.RSP_OK) {
            when (type) {
                1 -> dfuHandle.init()
                2 -> dfuHandle.sendPacket()
                3 -> dfuHandle.check()
                4 -> {
                    //The upgrade is successful, wait for the
device to restart
                    dfuHandle.endAndRelease()
                }
            }
        } else {
            //Upgrade exception or failure
        }
    }

    override fun onProgress(percent: Int) {
        // file upgrade progress
    }
}

```

## 2.3.7Manual measurement

```

//StarthHeartRateRsp parameter description
private byte type; type 1: heart rate 2: blood
pressure 3: blood oxygen
private byte errCode; measurement error code 0:
normal 1: measurement failed 2: measurement failed
private byte value; measurement value: heart rate
or blood oxygen
private byte sbp; blood pressure sbp

```

```

        private byte dbp; blood pressure dbp

        // manual heart rate
        BleOperateManager.getInstance().manualModeHeart(new
        ICommandResponse<StartHeartRateRsp>() {
            @Override
            public void onDataResponse(StartHeartRateRsp
            resultEntity) {

                }

            });
        //manual blood pressure
        BleOperateManager.getInstance().manualModeBP(new
        ICommandResponse<StartHeartRateRsp>() {
            @Override
            public void onDataResponse(StartHeartRateRsp
            resultEntity) {

                }

            });
        //manual blood oxygen
        BleOperateManager.getInstance().manualModeSpO2(new
        ICommandResponse<StartHeartRateRsp>() {
            @Override
            public void onDataResponse(StartHeartRateRsp
            resultEntity) {

                }

            });
        //manual pressure
        BleOperateManager.getInstance().manualModePressure(new
        ICommandResponse<StartHeartRateRsp>() {
            @Override
            public void onDataResponse(StartHeartRateRsp
            resultEntity) {

                }

            });
    };

```

## 2.3.9 Changes in Ring measurement data are proactively reported to the APP

```

//Add listener

BleOperateManager.getInstance().addOutDeviceListener(ListenerKey.He

```

```
art,myDeviceNotifyListener)
```

ListenerKey Parameter Description

```
public class ListenerKey {  
    public static int Heart=1;           Heart  
    public static int BloodPressure=2;   Blood Pressure  
    public static int BloodOxygen=3;     Blood Oxygen  
    public static int Temperature=5;     Temperature  
    public static int SportRecord=7;     Sport Record  
    public static int All=7;             All  
}
```

Monitoring instructions

```
    inner class MyDeviceNotifyListener : DeviceNotifyListener() {  
        override fun onDataResponse(resultEntity: DeviceNotifyRsp?)  
    {  
        if (resultEntity!!.status == BaseRspCmd.RESULT_OK) {
```

```
        BleOperateManager.getInstance().removeOthersListener()  
        when (resultEntity.dataType) {  
            1 -> {  
                //Watch heart rate test changes  
            }  
            2 -> {  
                //Watch blood pressure test changes  
            }  
            3 -> {  
                //Watch blood oxygen test changes  
            }  
            4 -> {  
                //Changes in watch step counting details  
            }  
            5 -> {  
                //Watch body temperature changes on the day  
            }  
            7 -> {  
                //Generate new exercise records  
            }  
            0x12 -> {  
                //7312 00005200025100003c0000000066  
                AwLog.i(Author.HeZhiYuan,  
ByteUtil.bytesToString(resultEntity.loadData))  
  
                val step = BLEDataFormatUtils.bytes2Int(  
                    byteArrayOf(  
                        resultEntity.loadData[1],  
                        resultEntity.loadData[2],  
                        resultEntity.loadData[3]
```

```

        )
    )

    val calorie = BLEDataFormatUtils.bytes2Int(
        byteArrayOf(
            resultEntity.loadData[4],
            resultEntity.loadData[5],
            resultEntity.loadData[6]
        )
    )

    val distance =
BLEDataFormatUtils.bytes2Int(
        byteArrayOf(
            resultEntity.loadData[7],
            resultEntity.loadData[8],
            resultEntity.loadData[9]
        )
    )

    deviceNotification(step, distance, calorie)
}
}
}
}
}

//Remove the listener. It must be removed after registration,
otherwise multiple callbacks will appear.

BleOperateManager.getInstance().removeNotifyListener(ListenerKey.He
art)
//Remove all listeners

BleOperateManager.getInstance().removeNotifyListener(ListenerKey.Al
l)

```

## 2.3.10 APP opens exercise type

```

// status 1 Start movement 2 Pause 3 Continue 4 End 6 Movement
start timestamp
//Sport type 4 Walking 5 Jumping rope 7 Running 8 Hiking 9
Cycling 10 Other sports 20 Hiking 21 Badminton
22 Yoga 23 Aerobics 24 Spinning 25 Kayaking 26 Elliptical
machine 27 Rowing machine 28 Table tennis 29 Tennis
30 Golf 31 Basketball 32 Football 33 Volleyball 34 Rock
climbing 35 Dance 36 Roller skating 60 Outdoor hiking
CommandHandle.getInstance().executeReqCmd(

```

```

        PhoneSportReq.getSportStatus(

            1, sportType
        ), gpsResponse
    )

private var gpsResponse: ICommandResponse<AppSportRsp> =
    ICommandResponse<AppSportRsp> { resultEntity ->
        AwLog.i(Author.HeZhiYuan, resultEntity)
        if (resultEntity != null) {
            when (resultEntity.gpsStatus) {
                6 -> {
                    //Exercise start time (Unit second)
                }

                2 -> {
                    //Exercise pause
                }

                3 -> {
                    // //Exercise continues
                }

                4 -> {
                    //Exercise end
                }
            }
        }
    }

//Report data during exercise
//Add motion data reporting and monitoring
BleOperateManager.getInstance().addSportDeviceListener(0x78,
myDeviceSportNotifyListener)
//Remove sports data reporting monitoring

BleOperateManager.getInstance().removeSportDeviceListener(0x78)
//Listening to inner classes
inner class MyDeviceNotifyListener :
DeviceSportNotifyListener() {
    override fun onDataResponse(resultEntity: DeviceNotifyRsp?)
    {
        super.onDataResponse(resultEntity)

        if (resultEntity!!.status ==
BaseRspCmd.RESULT_OK) {
            //Movement duration, unit seconds
            val sportTime = bytes2Int(

```



```

        val sportTime = bytes2Int(
            byteArrayOf(
                resultEntity.loadData[2],
                resultEntity.loadData[3]
            )
        )
    //Exercise real-time heart rate
    val heart = bytes2Int(
        byteArrayOf(
            resultEntity.loadData[4]
        )
    )

    //The number of steps generated during exercise
    will only have data when the exercise type is 4, 7, or 8, otherwise
    it will be 0
    val step = bytes2Int(
        byteArrayOf(
            resultEntity.loadData[5],
            resultEntity.loadData[6],
            resultEntity.loadData[7]
        )
    )

    //The distance generated during exercise will only
    have data when the exercise type is 4, 7, or 8, and the others will
    be 0.
    val distance = bytes2Int(
        byteArrayOf(
            resultEntity.loadData[8],
            resultEntity.loadData[9],
            resultEntity.loadData[10]
        )
    )

    //Calories generated during exercise
    val calorie = bytes2Int(
        byteArrayOf(
            resultEntity.loadData[11],
            resultEntity.loadData[12],
            resultEntity.loadData[13]
        )
    )

    //Error status during exercise
    val status = bytes2Int(
        byteArrayOf(
            resultEntity.loadData[1]
        )
    )
}

```

```

        ,
        val sportType = BLEDataFormatUtils.bytes2Int(
            byteArrayOf(
                resultEntity.loadData[0]
            )
        )
        if (status == 0x03) {
            //Not wearing
        }
    }
}

/**
 * Convert the byte array to int type, with the high byte of
the array first
 *
 * @param data
 * @return
 */
public static int bytes2Int(byte[] data) {
    int length = data.length;
    int res = 0;
    for (int i = 0; i < length; i++) {
        res |= (data[i] & 0xFF) << (8 * (length - 1 - i));
    }
    return res;
}

```